VOCO Grandio Core Dual Cure

(GB) **in accordance with DIN EN ISO 4049** **(CE)** ⁰⁴⁸²

Instructions for use

Product description:

Grandio Core Dual Cure is a dual-curing, radiopaque, flowable composite for core build-ups, with very good mechanical properties.

Grandio Core Dual Cure contains 77% w/w inorganic fillers in a methacrylate matrix.

Shades:

Grandio Core Dual Cure is available in three different shades:

- dentine shade for highly aesthetic demands, e.g., under all-ceramic restorations or highly translucent composite restorations;
- blue for reproducing the transition between the material and the dental tissue, as a base for metal crowns or faced crowns with a metal base
- white as an option for reproducing preparation borders, whilst simultaneously maintaining aesthetic requirements

Indications:

- Core build-up of vital and non-vital teeth

- Luting of root posts

Type of application:

Moisture control:

- Ensure that the working field is sufficiently dry. The use of a rubber dam is recommended.

Cavity preparation:

In principle, the cavities should be prepared according to the rules for adhesive restorative treatment, i.e. preparation should be minimally invasive in order to protect healthy dental hard tissue. Then clean and dry the cavity.

Lining the cavity:

Take appropriate measures in areas close to the pulp.

Bonding materials:

Grandio Core Dual Cure is applied with a dentine/enamel bond indicated for dual-curing composites, using the adhesive technique. We recommend the use of **Futurabond DC**. Please follow the respective instructions for use in regard to preparation and application.

Preparing the QuickMix syringe:

For technical reasons, there may be slight differences in the fill levels of the two syringe chambers prior to initial use. You should therefore remove the syringe cap and continue to press out the material until the material flows equally from both outlets. Then attach a type 11 mixing tip. The notches on the syringe and the mixing tip must correspond. The mixing tip is fixed by rotating clockwise through 90°. If necessary, an additional type 4 intrarotal tip can be placed on the end of the mixing tip. The material is automatically mixed in the correct proportions inside the mixing tip as it is dispensed and can be applied immediately. All QuickMix syringes are equipped with a removable plunger. The plunger can easily be pulled backwards out of the syringe and then reinserted. The seals in the syringe do not change their position whilst this is being done. This prevents contamination of the materials in the syringe.

Application of core build-up materials:

Place matrix around the prepared tooth or apply the shaping aid (**Rebildta Form**) and **Grandio Core Dual Cure** directly from the QuickMix syringe. **Grandio Core Dual Cure** has a working time of approx. 1.5 min (23°C).

Light-curing:

Apply mixed materials into the cavity in increments of no more than 2 mm and polymerise for 40 s each time using a LED/halogen lamp with a minimum light output of 500 mW/cm². Hold the light emitting window of the polymerisation lamp as close to the surface of the core build-up as possible, as the curing may otherwise be insufficient. The restoration can be finished immediately after light-curing.

Chemical curing:

If light-curing is not carried out, the build-up can be processed 5 min after insertion into the mouth at the earliest.

Luting of root posts:

Use an adhesive indicated for dual-curing composites for the luting of root posts. We recommend the use of **Futurabond DC**.

Please follow the respective instructions for use in regard to preparation and application.

Luting of root posts (in combination with Futurabond DC):

- For the luting of root posts, the mixed **Futurabond DC** (liquid 1 and 2) is massaged into the root canal and onto the occlusal surface of the core build-up for 20 s using an applicator. Dry the solvent with an air stream of oil-free air for 5 s. Remove any possible excess material in the canal with absorbent paper strips. Do not light-cure the bonding layer!
- It is recommended that **Grandio Core Dual Cure** is applied with the fine intrarotal tip type 1 directly into the root canal. Do not use Lentulo! In addition, **Grand Core Dual Cure** can be applied to wet the root post, when the latter has been prepared according to the manufacturer’s instructions. Insert the root post in such a way that slight amounts of excess cement are formed. Due to the increased temperature and the contact with the bonding agent, the working time of **Grandio Core Dual Cure** in the root canal is reduced to 40 s.
- To fix the root post and to cure the occlusally applied **Futurabond DC**, the material should be light-cured for at least 40 s.
- For core build-ups, proceed as described in "Application of core build-up materials".

Information/precautionary measures:

- Grandio Core Dual Cure** contains methacrylates, amines, benzoyl peroxide and BHT. Do not use in patients with known allergies to the constituents of these products.
- Phenolic substances, especially preparations containing eugenol and thymol, lead to curing failures in composites. The use of zinc oxide eugenol cements or other eugenol-containing materials in combination with composites is therefore to be avoided.
- The use of one-bottle adhesive systems is to be avoided, as an optimal bond with chemically or dual-cured materials cannot, as a general rule, be assured.
- Isolate the core build-up for taking impressions or fabricating temporary crowns/bridges (it may be sufficient to moisten it thoroughly), in order prevent any bonding with these materials.
- Avoid contact with the oral mucosa.
- In purely chemical curing, low temperatures (e.g., when working on the anterior region) can lead to an increase in setting time.

Storage instructions and application method:

After use, leave the mixing tip on the QuickMix syringe or replace the cap. Do not expose **Grandio Core Dual Cure** to direct sunlight or the operating light. Store at temperatures of 4°C to 23°C. Storage in a refrigerator is recommended. Bring to room temperature before use. Do not use after the expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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